

ANIOSYME DD1

Cleaning and
pre-disinfection
of instrumentation



- Patented formula
- First tri-enzyme liquid detergent with disinfecting properties
- Especially developed for collecting of soiled instrumentation: no crossed contamination

INDICATIONS

- Reinforced cleaning and pre-disinfection of medico-surgical instrumentation, medical devices and endoscopic equipment.
- Cleaning in ultrasonic bins.
- Collecting of soiled instrumentation.


CHARACTERISTICS

- First enzymatic liquid detergent with disinfecting properties.
- Tri-enzyme complex : protease, lipase, amylase, combined with surfactants
- Proved efficacy (DCP*) *Degreasing Cleaning Power
- Chlorine free formulation: no oxydization of materials.
- Aldehyde free formulation: no fixation of proteins.
- Stability of soaking bath: 8 hours (operating working day).
- Enzymatic stability proved.
- Neutral pH: compatible with alloys.
- Absence of chloride (ammonium propionate): non corrosive effect with materials.
- Compatible with ultra-sonic process.
- Clear liquid, blue perfumed

PRODUCT
TO BE DILUTED
25 ml > 5L




0,5%
5 min.



ANIOSYME DD1

Cleaning and pre-disinfection of instrumentation


INSTRUCTIONS FOR USE



1


0.5 % dilution:

Pour a 25 ml dose in 5 litres of cold or tepid water.




2

Totally immerse the medical device. **Advised soaking time: 5 minutes.** Brush when necessary. For endoscopic equipment: brush channels.



3

Rinse thoroughly with tap water (of good microbiological properties) For endoscopic equipment: both internal and external parts of the medical device.



4

Dry with a clean towel. Proceed to the next step (see protocol established by the service).

COMPOSITION

Quaternary ammonium propionate, polyhexamethylene biguanide hydrochloride, enzymatic complex (protease, lipase and amylase), surface-active agents, stabilising agents, sequestering agents, fragrance, colouring.

PRECAUTIONS FOR USE

Dangerous – respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: from +5°C to +35°C. Class IIb medical device (Directive 93/2/EEC as amended).

PACKAGING

- 1 Box of 200 dosis, 25 ml each..... Ref. 1200.097
- 2 12 dosing bottles, 1 litre each..... Ref. 1200.095
- 3 4 cans, 5 litres each with a 25 ml dosing pump..... Ref. 1200.036

MICROBIOLOGICAL PROPERTIES

Active against	Standards	Contact time
Bacteria	EN 1040, EN 13727 (dirty conditions : albumine and erythrocytes MRSA [EN 13727])	5 minutes
Mycobacteria	Mycobacterium tuberculosis (Tb)	5 minutes
Yeasts	EN 1275, EN 13624	5 minutes
Viruses	HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Herpes virus, Influenza virus A [H1N1], Vaccinia virus	5 minutes

